

HAFFKINE BIO PHARMACEUTICAL CORPORATION LIMITED

Procurement Cell

(A Government of Maharashtra Undertaking)

Regd. Office: Acharya Donde Marg, Parel, Mumbai 400 012 (INDIA)

Phone No: 022- 24129320-23

Managing Director: 022-24150628

General Manager (Procurement Cell):022-

24100478

Website: http://www.vaccinehaffkine.com

E-mail: procurementcell@vaccinehaffkine.com

No.: 2364 / Haffkine / Procurement Cell/E-222/C-Arm

Compatible Table/2019-20

Date: 26/12/ 2019

To.

M/s. Staan Bio-Med Engineering Pvt. Ltd.,

190-A, Bharathiar Road, Ganapathy,

Coimbatore - 641006.

E-Mail: Stanly@staan.in

Subject: Supply Order for Tender No. E-222/ C-Arm Compatible Table

Reference: 1. Tender No. E-222/HBPCL/PC/ C-Arm Compatible Table /2017-18

2. Sanction of Tender Approval Committee Meeting Dated 10.12.2019

With reference to the tender cited under reference no 1, you are requested to supply the following goods as per details mentioned below to consignee list enclosed with this order.

Sr.No	Name of the item	Specification of item	Quantity/Unit (DMER)	Unit Rate inclusive of GST(Rs.)	Total Amount Rs.
1	C-Arm Compatible Table Model : Staan & Glory 12 Function	As per Annexure X	2	9,73,000/-	19,46,000/-

Total amount in words: Rupees Nineteen Lakhs Forty Six Thousand Only.

Factory Location: M/s. Staan Bio-Med Engineering Pvt. Ltd., 190-A, Bharathiar Road, Ganapathy, Coimbatore – 641006.

- **Forwarding:** Forwarding Free on Road Destination. I.e. door delivery basis The delivery should be damage free & intact Packaging.
- Delivery Period: 06 weeks from the date of receipt of order by the supplier to the consignee attached.
- Pre-Dispatch Inspection: Supplier shall make necessary arrangement / facilitate to carry out Pre-Dispatch inspection as per Tender Terms & condition and submit the Inspection report to this office. The Pre-Dispatch inspection cost will be borne by supplier. Machine should be dispatched only after satisfactory Pre-Dispatch Inspection.

- Risk purchase clause: If the bidder fails to supply the stores within the stipulated delivery period, the order will stand cancelled. Undersigned shall be entitled to purchase such stores from any other source at such price which ordinarily should not be more than 10% of the tender price. The extra expenditure in such cases shall be recovered by Managing Director, Haffkine Bio Pharmaceutical Corporation Ltd. (Procurement Cell), Mumbai from the Supplier.
- Payment Terms: Initial Payment of 90% of the contract value will be made within 8 weeks on delivery and successful installation. Remaining 10% of the contract value will be paid only on satisfactory commissioning and operating of the machinery.
- 6 Labelling: The word "For use of GOVERNMENT OF MAHARASHTRA NOT FOR SALE" should be printed on each unit pack in readable Purple or Green Colours.
- 7 Acceptance & Receipt: It should be submitted in Appropriate Format to the purchasing authority.
- 8 **Delivery Challan** Should be sent in the name of consignee in duplicate. It should specify Name of Equipment / Mfg. by / packing & quantity.
- 9 Invoice Copy Should be sent in triplicate on the Name of Managing Director, Haffkine Bio Pharmaceutical Corporation Ltd.(Procurement Cell), Mumbai. Along with Bill of Entry and Country of Origin Certificate of the consignment.

10 Other Terms:

- 1) Warranty: The warranty period shall be for 2 years from the date of commissioning of all equipment supplied as certified by the consignee. After completion of 2 years warranty period Manufacturer/Supplier should give commitment to ensure services and supply of spare part for further 8 years. The successful tenderer must ensure 100% uptime during warranty period. In case of downtime, warranty period will be extended for period of downtime. If the equipment is not attended within 24 hours for Mumbai and 48 hours for other places the supplier will be liable to pay a penalty of 0.07% of purchase cost for every day of delay. Such penalty will be recovered from the amount of security deposit. Certificate of such uptime / downtime issued by the end user will be binding for the supplier Replacement of spares parts thereof due to manufacturing defects during warranty period will be entirely at the supplier's cost.
- 2) The user institution will enter to the Comprehensive Maintenance Contract with supplier agency @ 5% of the Written down value of the equipment per year for 8 years after completion of warranty period. In case of non-compliance of CMC the supplier will be liable to pay penalty or for appropriate action. Payment of CMC on yearly basis will be made by the user's institution, at the end of the year after satisfactory performance report from the end user.
- Contract Agreement: Bidder should submit Contract Agreement on non-judicial stamp paper of requisite value.

Fall Clause

It is a condition of the contract that all through the currency thereof, the price at which you will the supply stores should not exceed the lowest price charged by you to any customer during the currency of the contract and that in the event of the prices going down below the rate contract prices you shall promptly furnish such information to us to enable to amend the contract rates for subsequent supplies.

12

The Bidder should submit (within 7 days) amount of 1.5% ie. Rs. 29,190/- of order value to meet other incidental expenditure and 3% i-e Rs. 58,380/- as Security Deposit in form of Bank Guarantee. The Bank Guarantee Should be Valid for 8 Years after the expiry date of warranty, issued by any Nationalized/Scheduled Bank

Amount to be deposited to Following Account:

Name of Account	Haffkine Bio-Pharmaceutical Corporation Ltd. (Procurement Cell), CESS Account Mumbai.	
Name of the Bank & Branch	Canara Bank, Branch-Parel	
Account No.	0110201004893	
IFSC Code	CNRB0000110	

Consignee: As per list enclosed

मा. व्यवस्थापकीय संपालक यांच्या मान्य टिपणी नुसार

ruth.

(Dr. Rajesh Deshmukh)
Managing Director
Haffkine Bio Pharmaceutical Corporation Ltd.
(Procurement Cell), Mumbai.

General Manager (Procurement) HBCL Mumbai

Copy to:

- 1) Commissioner Health Services, Mumbai.
- 2) Director, Medical Education & Research, Mumbai-400 001.
- 3) Account Manager, Haffkine Bio Pharmaceutical Corporation Ltd. (Procurement Cell), Mumbai.
- 4) Office File.

<u>Copy to</u> Consignee: Government Medical College & Super Specially Hospital, Nagpur. : As per Tender Condition No.17 The user Institution should get the Comprehensive Maintenance Contract done with supplier agency @ 5% of the Written down value of equipment per year for Eight years after Completion of warranty period.

Copy Submitted to: 1) Secretary, Medical Education & Drug Department, Mantralaya, Mumbai.

Annexure-X

	Technical specifications of C-Arm Compatible Table
A	OT Table should be a C-Arm Compatible Table, non hydraulic leak proof design, Fully Automatic Remote Controlled with International Quality Certifications – US FDA/ European Directive CE Certificate from European Notified Body (Certificate to be submitted)
>	It should be manufactured in ISO 9001 and ISO 13485 Certificate Facility
>	Table should be very sturdy with Load capacity – 150 kg
>	Design:
	Table should be made of SS 304 Construction and imported integrated hard plastics – Water % Spill Proof for long life and rust – free
>	The table should be single pillar eccentric design where the weight should come at the middle section to avoid cantilever issues.
>	The table base must have a non-metallic mounded plastic exterior for easy maintenance & for dent free long life.
>	All lower limb and upper body crisp images should be easily obtained with C-Arm with good stability.
A	Table top should have4 section radiolucent table top with very low density for crisp C-Am images.
^	This should be made of a composite material, crevices free, with no Screws/ holes on the table top surface, to prevent infection due to blood getting accumulated in the crevices, which is difficult to clean.
~	Side Rails: sturdy Side Rails fitted on both sides for attachments
>	Dimensions:
	Width: Table top should have a normal width of 480mm. For Future upgradeability, The table top should have provision for detachable width extenders on both sides, complete with paddin and side rails, to achieve extra width of minimum 100 mm
	Controls:
>	It should have dual control console – One on the table and the other via corded remote control Remote control to have clear indications for all different positions
^	Cordless Remote control: The table must also have an independent cordless remote control handset operating on Bluetooth technology, so that there is no interference to other equipment inside the operation theater. This should facilitate the anesthetist to operate the table from a distance without causing interference to the surgeon and his team.
>	Emergency override control: the table must have an independent control system apart from control handsets so that it should be possible to operate the table even without the remote handset unit.

		1 echnical specification	ons of C-Arm Compatible Table			
	All ele	ctrical operations and moveme	ents to be very smooth, jerk free using top quality			
	electrical actuators.					
_	Auto level Functions: Tilts of the table can be brought back to neutral level flat zero position					
>	with press of button and remote control handset.					
	Constr	uction:				
>	No hol	es, nut bolts on table top thus	avoiding accumulation of blood and stains			
>	Head rest: Should be Detachable					
>	Clamps & Brackets: All clamps for brackets and hand locks must be made of stainless steel					
	Mobili	ty: fitted with castor wheels for	or free table mobility & floor lock. The floor lock for			
>	castor	wheels must have a central loo	cking facility with a single operation.			
>	Kidney Bridge should be detachable and not fixed – so that there is no impediment to C-Arm imaging on the table top and can be fixed anywhere on the table top.					
	arrange	ement, so that leg position can	we normal lithotomy position with a brake & lock be easily maneuvered single handedly without it falling			
			d out horizontally, giving you wide & easy access between a down as desired or extended out as per requirement.)			
>	the leg	s. Legs can also be tilted up &				
^	the leg	s. Legs can also be tilted up & should be upgradeable with all TECHNICAL	down as desired or extended out as per requirement.)			
	the leg	s. Legs can also be tilted up & should be upgradeable with all TECHNICAL SPECIFICATIONS	Il future attachments so that it will not get obsolete.			
	Table s	s. Legs can also be tilted up & should be upgradeable with all TECHNICAL SPECIFICATIONS Table Top Dimensions	2 down as desired or extended out as per requirement.) 21 future attachments so that it will not get obsolete. 2150 mm Length × 480 mm Width			
	Table s	s. Legs can also be tilted up & should be upgradeable with all TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections	2 down as desired or extended out as per requirement.) 21 future attachments so that it will not get obsolete. 2150 mm Length × 480 mm Width 4 nos (including Headrest)			
	Table s	s. Legs can also be tilted up & Should be upgradeable with all TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment	2 down as desired or extended out as per requirement.) 2 If future attachments so that it will not get obsolete. 2 150 mm Length × 480 mm Width 4 nos (including Headrest) 8 40 mm to 109mm			
	Table s 1 2 3 4	TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide	2150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm			
	Table s	s. Legs can also be tilted up & Should be upgradeable with all TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg	2150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm 30 Degrees			
	Table s 1 2 3 4 5	TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide	2150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm			
	1 2 3 4 5 6	TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg	2 down as desired or extended out as per requirement.) 2 If future attachments so that it will not get obsolete. 2 150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm 30 Degrees 30 Degrees			
	1 2 3 4 5 6 7	TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg Lateral Tilt	2150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm 200 mm 30 Degrees 30 Degrees 20 Degrees			
	1 2 3 4 5 6 7 8	s. Legs can also be tilted up & Should be upgradeable with all TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg Lateral Tilt Flex and Reflex	2 down as desired or extended out as per requirement.) 2 150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm 30 Degrees 30 Degrees Flex is 80 Degrees and Reflex is 45 Degrees By press of a single button on both consoles to bring back zero position from all tilts			
	1 2 3 4 5 6 7 8	s. Legs can also be tilted up & Should be upgradeable with all TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg Lateral Tilt Flex and Reflex	2 down as desired or extended out as per requirement.) 2 I future attachments so that it will not get obsolete. 2 150 mm Length × 480 mm Width 4 nos (including Headrest) 8 40 mm to 109mm 2 00 mm 3 0 Degrees 3 0 Degrees 2 0 Degrees Flex is 80 Degrees and Reflex is 45 Degrees By press of a single button on both consoles to bring back			
	1 2 3 4 5 6 7 8 9	TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg Lateral Tilt Flex and Reflex Auto Level Feature	2 down as desired or extended out as per requirement.) 2 150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm 30 Degrees 30 Degrees Flex is 80 Degrees and Reflex is 45 Degrees By press of a single button on both consoles to bring back zero position from all tilts			
	1 2 3 4 5 6 7 8 9 10	TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg Lateral Tilt Flex and Reflex Auto Level Feature Emergency Override	2 down as desired or extended out as per requirement.) 2 I future attachments so that it will not get obsolete. 2 150 mm Length × 480 mm Width 4 nos (including Headrest) 8 40 mm to 109mm 2 00 mm 3 0 Degrees 3 0 Degrees Flex is 80 Degrees and Reflex is 45 Degrees By press of a single button on both consoles to bring back zero position from all tilts Emergency Override feature for Remote control			
	1 2 3 4 5 6 7 8 9 10 11	TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg Lateral Tilt Flex and Reflex Auto Level Feature Emergency Override	2150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm 30 Degrees 30 Degrees Flex is 80 Degrees and Reflex is 45 Degrees By press of a single button on both consoles to bring back zero position from all tilts Emergency Override feature for Remote control Table on 230 V AC.			
>	1 2 3 4 5 6 7 8 9 10 11	TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg Lateral Tilt Flex and Reflex Auto Level Feature Emergency Override Supply Voltage	2 down as desired or extended out as per requirement.) 2 150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm 30 Degrees 20 Degrees Flex is 80 Degrees and Reflex is 45 Degrees By press of a single button on both consoles to bring back zero position from all tilts Emergency Override feature for Remote control Table on 230 V AC. Motors should operate on 24v DC			
>	1 2 3 4 5 6 7 8 9 10 11	s. Legs can also be tilted up & Should be upgradeable with all TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg Lateral Tilt Flex and Reflex Auto Level Feature Emergency Override Supply Voltage rd Accessories: Detachable flat Head Plate	2 down as desired or extended out as per requirement.) 2 150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm 30 Degrees 20 Degrees Flex is 80 Degrees and Reflex is 45 Degrees By press of a single button on both consoles to bring back zero position from all tilts Emergency Override feature for Remote control Table on 230 V AC. Motors should operate on 24v DC			
>	1 2 3 4 5 6 7 8 9 10 11 Standa	TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg Lateral Tilt Flex and Reflex Auto Level Feature Emergency Override Supply Voltage rd Accessories:	2 down as desired or extended out as per requirement.) 2 150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm 30 Degrees 20 Degrees Flex is 80 Degrees and Reflex is 45 Degrees By press of a single button on both consoles to bring back zero position from all tilts Emergency Override feature for Remote control Table on 230 V AC. Motors should operate on 24v DC			
>	1 2 3 4 5 6 7 8 9 10 11 Standa	TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg Lateral Tilt Flex and Reflex Auto Level Feature Emergency Override Supply Voltage rd Accessories: Detachable flat Head Plate Adjustable Body Side Supp	2150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm 30 Degrees 20 Degrees Flex is 80 Degrees and Reflex is 45 Degrees By press of a single button on both consoles to bring back zero position from all tilts Emergency Override feature for Remote control Table on 230 V AC. Motors should operate on 24v DC			
>	1 2 3 4 5 6 7 8 9 10 11 Standa	Should be upgradeable with all should be upgradeable with all TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg Lateral Tilt Flex and Reflex Auto Level Feature Emergency Override Supply Voltage rd Accessories: Detachable flat Head Plate Adjustable Body Side Supply Swivelling Arm Rests – 2nd Mattress 40 mm Thick – 1nd Detachable Kidney Bridge	2150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm 30 Degrees 20 Degrees Flex is 80 Degrees and Reflex is 45 Degrees By press of a single button on both consoles to bring back zero position from all tilts Emergency Override feature for Remote control Table on 230 V AC. Motors should operate on 24v DC			
>	1 2 3 4 5 6 7 8 9 10 11 Standa	Should be upgradeable with all should be upgradeable with all TECHNICAL SPECIFICATIONS Table Top Dimensions No of Sections Height Adjustment Table top Slide Trendelenberg Reserve Trendelenberg Lateral Tilt Flex and Reflex Auto Level Feature Emergency Override Supply Voltage rd Accessories: Detachable flat Head Plate Adjustable Body Side Supply Swivelling Arm Rests – 2nd Mattress 40 mm Thick – 1nd Detachable Kidney Bridge	2150 mm Length × 480 mm Width 4 nos (including Headrest) 840 mm to 109mm 200 mm 30 Degrees 20 Degrees Flex is 80 Degrees and Reflex is 45 Degrees By press of a single button on both consoles to bring back zero position from all tilts Emergency Override feature for Remote control Table on 230 V AC. Motors should operate on 24v DC			

Consignee Details

	M/s. Staan Bio-Med Engineering Pvt. Ltd.		
	C-Arm Compatible Table		
	Model: Staan & Glory 12 Function		
Delivery Period	06 weeks		
PO Reference No.	No.: 2364 /Haffkine/Procurement Cell/E-222/ C-Arm Compatible Table/2019-20 Date: 26 12 20 9		
Sr. No.	Name & Address of the Consignee		
1)	Government Medical College & Super Specialty Hospital, Medical Square Hanuman Nagar, Nagpur - 440003.	02	

मा. व्यवस्थापकीय संचालक यांच्या मान्य टिपणी नुसार

General Manager (Procurement) HBCL Mumbai (Dr. Rajesh Deshmukh)
Managing Director
Haffkine Bio Pharmaceutical Corporation Ltd.
(Procurement Cell), Mumbai.